


FLECOTEC

FLECOBOX (I-Box)



The table is optional and is available in various models for the I-Box

FLECOBOX I-Box

THIS INDIVIDUALLY PREFABRICATED FLEXIBLE ISOLATOR WITH SECONDARY CONTAINMENT AND INTEGRATED GLOVES MEETS ALL THE REQUIREMENTS OF CONTAINMENT LEVEL OEB 5 FOR HIGHLY TOXIC PRODUCTS.

The design of the I-Box does not require a base frame but is given its shape by the PE support tubes on the inside.

As a result and due to its very compact design, the I-Box is preferably used in laboratories where it offers personnel an additional high level of safety for a wide range of work steps.

STANDARD SIZES

	250	500	650	800
Width A (mm)	250	500	650	800
Length B (mm)	450	500	650	650
Height C (mm)	450	500	650	650
Base outlet possible	–	✓	✓	✓
Number of gloves (n)	2	2-4	2-4	2-6

(special sizes on request)



Flecotec AG

Breisgauring 13
79427 Eschbach
Germany

T + 49 7632 828 553-0
F + 49 7632 828 553-9

mail.ft@rommelag.com
www.rommelag.com

